


SUPiA gram Viewer

Able to check image **anywhere** you want



Signers CR or DR



SUPiA gram
Stellar gram











▶ Easy to access patients' images through Signers CR or DR



Make ID and PW in SUPIA gram or Stellar gram
And Set IP address of SUPIA gram or Stellar gram PC

SUPiA gram Viewer

▶ Able to write patients' images with simple viewer to CD or USB

 D3DCompiler_43.dll	2020-04-21 오후 6:54
 d3dx11_43.dll	2020-04-21 오후 6:54
 Dicom32.dll	2020-04-21 오후 6:54
 FillPixelShader.hlsl	2020-04-21 오후 6:54
 FillVertexShader.hlsl	2020-04-21 오후 6:54
 TexPixelShader.hlsl	2020-04-21 오후 6:54
 TexVertexShader.hlsl	2020-04-21 오후 6:54
 Viewer	2020-04-21 오후 6:54

Execute 'Viewer' file in app folder
Images will pop up automatically like right pic.



▶ Simple Tools



- ✓ Measurement of length
- ✓ Angle
- ✓ Squaremeter
- ✓ Annotation
- ✓ Cobbs angle
- ✓ Text
- ✓ Window W/L
- ✓ Back to original

SUPiA gram Viewer

► Powerful Print Option

- 1) Adding up Hspital Logo (PNG)
- 2) Rotating Images
- 3) W/L
- 4) Magnifying Images

