

# CNME011213


12-channel Electrocardiograph



## Features

- 5.7 inch TFT color screen
- Lead status indication
- Simple & clear keyboard layout, back light support
- Internal memory of 300 ECG records
- High accurate digital filter (EMG filter, AC filter, Drift filter, lowpass filter)
- 122 kinds of diagnosis report
- 120 seconds ECG waveforms review and print
- Support 210mm rolling, Z-folded paper
- USB port for external printer connection & data transmission
- PC-ECG management system
- Working mode: Auto, Manual, Rhythm

## Interface



## Various printing mode



Printing mode: 12X1, 6X2+1R, 3X4+1R, auto analysis



# CNME011213

12-channel Electrocardiograph



**NANJING KAIHONG HEALTHCARE CO., LTD**

ADD: No.101 Caochangmen Street, Gulou District, Nanjing, China.  
TEL: +86(0)25 83617669 +86(0)13611512069  
Email: meditech@cn-meditech.com elvanenglish1984@gmail.com  
WEB: www.cnmeditech.com www.kh-medical.com





# CNME011213

12-channel Electrocardiograph



Hidden handle




Alphanumeric keyboard;  
Back light for button



Memory of 300 ECG records;  
USB for data transmission;  
PC-ECG management system

## Clear professional report



12X1 + Average template + Automatic analysis

✓ 5.7 inch TFT color screen,  
easy for observation

✓ Internal aluminum alloy rack,  
shockproof, Anti-interference  
design

✓ High quality printing  
system, print is clear and  
smooth

✓ Internal modular design,  
easy for maintenance

