



ISO 13485



NORTHERN
MEDITEC

Committed to Excellence



NORTHERN

Northern Meditec Limited

ADD: 4th Floor, Building C, Jin Wei Yuan Industrial Park, Julongshan Area,
Kengzi St, Pingshan District, Shenzhen, P. R. China

TEL: +86 755 29970566 FAX: +86 755 86528647


EMAIL: info@northernmeditec.com WEBSITE: www.northernmeditec.com

PRODUCT CATALOGUE

About Northern



Northern meditec limited is a professional manufacturer of high-end medical equipment. Since established, Northern has been keeping the value of committed to excellence, to R&D and manufacture the best medical devices such as patient monitor, anesthesia machine, ventilator, infant radiant warmer incubator, etc.



Anesthesia Patient Monitor Gemini

Unique Design



Features

- 15.6"/19" switchable TFT LCD touch screen, wide and flat design
- Full-touch screen control, no physical buttons at front panel
- Aluminum material shell with effective heat dissipation capability
- Fanless design allows for quiet care environment
- Anti-motion and low perfusion SPO2 measurement
- Built-in rechargeable battery



Multiple-parameter options & Flexible screen size options

Critical Care Patient Monitor Venus



Features

- 15.6" TFT LCD Touch screen, wide and flat design
- ST and Arrhythmia analysis
- Big font, ECG 7-lead half-screen, ECG 7-lead full screen, ECG 12-lead full screen, OxyCRG, DynaTrend Multi-interface display
- 10 waveform display, up to 12-lead ECG analysis
- Powerful calculation (MEWS, Hemodynamic, Dose, Oxygenation, Ventilation, Renal function)
- Anti-motion and low perfusion SPO2 measurement

Up to 10 waveforms display



Optional plug-in module (plug and play)



Various external ports



15.6" TFT touch screen with big knob

PATIENT MONITORING

12.1" Multi-para Patient Monitor Taurus / Pisces / Virgo



Features

- 12.1" TFT LCD screen, wide and flat design
- Big font, ECG 7-lead half-screen, ECG 7-lead full screen, ECG 12-lead full screen, OxyCRG, DynaTrend Multi-interface display
- 8 waveform display, up to 12-lead ECG analysis
- Anti-motion and low perfusion SPO2 measurement
- Built-in rechargeable battery

Bull horn design handle



Optional plug-in module (plug and play)



Vital Sign Monitor Aquarius



Masimo SET® SpO2



Ear Infrared Thermometer



Bar/QR Code Scanner

Features

- 8.4" TFT color LCD, touch screen is optional
- Up to 6 parameters or 3 waveforms display
- Optional Ear/Forehead Infrared Thermometer, QR Code Scanner, Thermal Recorder
- Powerful calculation (MEWS, Hemodynamic, Dose, Oxygenation, Ventilation, Renal function)
- Anti-motion and low perfusion SPO2 measurement
- Built-in rechargeable battery (5 hours uninterruptable working time, can upgrade to 15 hours)



Portable Design




Touch Screen (Optional)



Emergency bag (Optional)

Central Monitoring System CMS



Features

- Networked data management of all Northern brand monitors
- Dual-screen display monitoring
- Connection via LAN/WLAN up to 64 bedside monitors
- Support WIFI/3G/4G connection, Cloud servers



Veterinary Monitor Taurus Vet / Aquarius Vet



Taurus Vet (12.1")



Aquarius Vet(8,4")

Features

- Suitable for feline, canine and other animals
- Multi-selectable user screen options
- Graphical & tabular trend review(120 hours)
- Built-in rechargeable battery
- Big font and color font display setting
- Spot-check and continuous monitoring mode
- Thermal recorder support



Fully electronic Anesthesia machine Atlas N7

Fully Electronic Control and Display



15.6" Touch screen




Auto FiO₂



Features

- 15.6" TFT touch screen, integrated patient monitor design(optional)
- Tidal volume:10ml-1500ml
- All ventilation modes to satisfy neonate, pediatric and adult patients
- Upto 4 waveforms display, 3 loops diagram
- User friendly design interface, easy to operate
- Double touch screen, special and simple-operatinon screen layout
- Human engineering design, rotatable supporting arm with adjustable angle.
- Full-electronic control and electronic flow meter, precise control and accurate data monitoring
- Autocalvable Metallic material of breathing system
- Integrated breathing circuit with built-in heater

Semi-electronic Anesthesia Machine Atlas N5



12.1" Touch screen



Plug-in module



AGSS



Auxiliary Outlet

Features

- 12.1" TFT touch screen
- All ventilation modes to satisfy neonate, pediatric and adult patient
- User friendly design interface
- Mechanical control mixer, electronic flowmeter
- Autocalvable Metallic material of breathing system
- Integrated breathing circuit with built-in heater

