

TECHNICAL DATA SHEET OPERATING TABLE ST1

Manufacturer	Sordina S.r.l.
Model	ST1
CND Code	Z12011202
RDM Code	1732518
Release on the market of the latest version (year)	2016
Intended Use	The Operating Table ST1 is designed to be used in operating rooms and for medical purposes, according to current legislation, to perform operations of all surgical specialties. The Operating Table ST1 is intended to be used in public and private healthcare environment and by personnel trained and authorized by the user (responsible authority). The Operating Table ST1 is suitable for adult patients.
Materials	
Base	Stainless Steel
Coloumn	Stainless Steel
Table	Stainless Steel
Lateral Rails for accessories	Stainless Steel (mis. 25 x 10 mm, standard in compliance with current European legislation). Rails are present for the whole lenght of the operating table.
Mattress	Mattresses in foamed polyurethane, radiolucent, free from latex, with antistatic treatment (in compliance with current European legislation). Mats covered in waterproof material and resistant to trauma, easily removable for disinfection and cleaning. Innovative ergonomic shape designed for antidecubitus.
Table Width	570 mm (external rails) - 520 mm (internal rails)
Table Lenght	2020 mm
Loading Capacity	297 Kg
Base features	Mobile type base on wheels, equipped with an innovative electrically operated floor locking system (electro-hydraulic locking).

Sordina S.r.l.

Capital € 10.000,00

P. IVA e C.F: 02978720593 - REA LT 214450

e-Mail: info@sordina.com-web: www.sordina.com

Administrative Office

Galleria del Pozzo Rosso nr.13 - 36100 Vicenza, Italy
Tel. +39 0444 233750 – Fax +39 0444 233790

Registered Office and Operating Office

Via dell'Industria 1/A - 04011 Aprilia (LT), Italy
Tel. +39 06 92730650

	Wheels made of materials and dimensions that ensure maximum mobility and safety, the absence of prints on the antistatic floors of the operating rooms and make the movement of the Operating Table easier.
Coupling sections	Coupling of the different sections of the operating table through an innovative lock-unlock system. The new system makes the assembly and disassembly of the various sections of the operating table extremely simple and fast.
Automatic electro-hydraulic movements of the operating table	Height Variation: 680-940 mm
Manual movements of the Operating Table	Height Variation: 680-940 mm Trendelemburg-Reverse: + 30 ° / -30 ° Lateral tilt (Tilt): + 20 ° / -20 ° Backrest Adjustment: + 75 ° / -25 ° Longitudinal translation: 300 mm Kidney Bridge: 136 mm Opening legs: 180 ° Leg Section Adjustment: 5 ° / 90 °
Easy disinfection and washing	Washing with traditional water spray systems. Waterproof and easily removable pillows for disinfection and cleaning.
CE Marking	Yes, according to 93/42 / CE directive The manufacturer is certified according to the UNI EN ISO 9001 and UNI CEI EN ISO 13485 standards for the design, production, marketing, technical assistance and installation of operating tables with accessories.

Sordina S.r.l.

Capital € 10.000,00

P. IVA e C.F: 02978720593 - REA LT 214450

e-Mail: info@sordina.com-web: www.sordina.com

Administrative Office

Galleria del Pozzo Rosso nr.13 - 36100 Vicenza, Italy
Tel. +39 0444 233750 – Fax +39 0444 233790

Registered Office and Operating Office

Via dell'Industria 1/A - 04011 Aprilia (LT), Italy
Tel. +39 06 92730650